overfibers
improve your performance

OverCEM SA

**Self-Adhesive
Innovative Formula**

Strong


Dual-curing

**Auto-mixing
T-Mixers tips**

**Multi-purpose
Luting cement**

**SIMPLE
STRONG
RELIABLE**

Made in Italy

		OverCEM SA		
Application		Universal	Translucent	Opaque
Self Adhesive		X	X	X
Shade		Yellow (A2)	Neutral	Dull White
Metal	Inlays / Onlays / Partial crowns	X		X
	Crowns / Bridges	X		X
	Ceramic/Metal Ti-Base Implant Abutments			XX
	Cast Posts and dowels	X	X	X
	Maryland Bridges	X		X
Glass-Ceramics	Inlays / Onlays / Partial Crowns	X	X	X*
	Veneers	X	X	X*
	Crowns	X	X	X*
Lithium-Disilicate	Inlays / Onlays / Partial Crowns	X	X	X*
	Crowns	X	X	X*
	3-units bridges	X	X	X*
	Veneers	X	X	X*
Oxide-Ceramics (Zirconia /Alumina)	Crowns, Bridges, Adhesive Bridges, Veneers	X	X	X*
	Ceramic/Zirconia Implant Abutments	XX	XX	X
	Endodontic Posts	X	XX	
Composite materials	Inlays / Onlays, Veneers	X	X	X*
	Crowns	X	X	
	Endodontic Glass-Fiber Posts	X	XX	

Tab. 1 X = recommended, XX = strongly recommended
X* = recommended in case of dark substrates

OverCEM SA

Dual-Curing Self-Adhesive Resin Cement

Instructions for use

Description of the Product

OverCEM SA is dual-curing, self-adhesive resin luting cement. It contains acidic adhesive monomers that create a strong bonding between the mineral part of dental tissues and the metal oxides of various restorative materials. The double syringe-type package is provided with small auto-mixing tip system (T-Mixer and T-Mixer endo) to reduce the cement waste inside tips.

OverCEM SA is available in three versions (Tab. 1), allowing a wide choice for any clinical or technical requirement:

OverCEM SA:

1. Universal, A2 colored;
2. Translucent, for fiber post cementation and neutral shading effects;
3. Opaque, when a masking effect of the substrate is necessary;

Purpose and fields of application of OverCEM SA

Strong adhesive connections of various restorative materials in implant and conventional prosthodontics, such as: metal/metal, metal/ceramic, metal/composite, metal/zirconia, and tooth substrates. OverCEM SA neither require dental tissue conditioning nor adhesive systems. OverCEM SA is distinguished by outstanding adhesive strength, different self-polymerizing times, and a reactive light-polymerization capability.

Composition

OverCEM SA contains mono and bi-functional (meth)acrylates. The percentage of inorganic filler is about 65% in weight and the dimension of the filler particles (D50) are about 5 µm. The

volumetric mixing ratio of base and catalyst pastes is 1:1.

Precautions

For the patient and for the dental office personnel:

- During the use of the cement both dental operator and patient must wear safety goggles to prevent any possible contact with eyes. In case of contact with eyes, wash immediately and abundantly with water and consult a doctor. For better results, use the mixing tips; it is not advisable to mix the cement after extruding the two pastes on the mixing plate. If the OverCEM SA cement is stored at low temperature, take it out at least 1 hour before using it.

Do not use OverCEM SA on substrates contaminated by temporary cements containing phenols (i.e. eugenol) since they may inhibit the polymerization.

- Do not use hydrogen peroxide (H₂O₂) because the free oxygen has a negative effect on both adhesion and polymerization of OverCEM SA.

For the patients:

- Avoid the use of this product on patients with known allergies to acrylates and peroxides.
- In case of prolonged contact with soft oral tissues, rinse abundantly with water and if it is possible with sodium bicarbonate in water solution. If an allergic reaction occurs, remove the product from the mouth as soon as possible, then consult a doctor.

For the dental office personnel:

- It is advisable to use protective gloves avoiding any skin contact of the product with uncured cement in order to minimize the risk of irritation and allergic reactions. In case of accidental contact with the skin, it is advisable to wash the part with soap and warm water.

To avoid damages to instrument, devices or surfaces, remove the cement before the polymerization takes place. Use acetone or ethanol to remove the uncured cement.

• If an allergic reaction occurs, consult a doctor. Material Safety Data Sheet of OverCEM SA are available upon request at www.overfibers.com, or contact your retailer.

Pulp protection

Not necessary if the distance from the pulp is >0.5mm. Provide pulp protection with fast-curing calcium hydroxide if the distance is less than 0.5mm.

Tooth preparation

Preparation of the cavity/tooth abutment

Before cementing

1. Clean carefully the prepared abutment or the dental surface with rotary brush and one solution composed of Clorexidina 0,2% + Cetrimide 0,2%. Wash thoroughly with air/water spray, then air-dry making sure the air is free of water or oil.
2. Do not dehydrate the dentin: an excessive drying can cause poor cement reaction and post-operative sensitivity.
3. Do not use substances such as EDTA, desensitizers, astringents, dentin sealants, or mouthwashes on the surfaces to be cemented. It is advisable to etch the enamel with phosphoric acid 32-37% for 15-20 seconds and remove it completely with water/air spray for at least 20 seconds.

WARNING! Do not acid-etch the dentin.

Post cementation: preparation of the root canal

1. Prepare the root canal with the state-of-the-art endodontic techniques.
2. Maintain an apical seal of at least 3-4 mm.
3. Rinse the root canal with a solution composed of Clorexidina 0,2% + Cetrimide 0,2%. Avoid any acid-etching inside the root canal. Avoid any hydrogen peroxide based disinfectants.
4. Rinse with water and dry with paper tips.
5. The use of the rubber dam is recommended.

WARNING! Do not use canal sealants containing eugenol as the polymerization of OverCEM SA may be hindered.

General recommendations for adhesive prosthodontics - Maryland bridges and laminate veneers

To guarantee a good adhesion, the abutment teeth must be sound or correctly restored, they should be in good periodontal conditions and have plenty of enamel surface available for cement adhesion. If dentin areas are exposed, they should not be acid etched. The prosthesis materials must be treated in accordance with the manufacturer's indications. The technique and the treatment used are under the responsibility of the dentist.

During cementation, apply the rubber dam to avoid the contamination of the surfaces. Etch the prepared enamel surfaces of the teeth for 15-20 seconds with 37% orthophosphoric acid; rinse thoroughly with water spray and dry with air free of water or oil.

Preparation of the restoration surfaces

OverCEM SA is designed to bind to the mineral component of the tooth. Any foreign substance placed between the cement and the dental structure hinders its adhesion. In particular, the use of acids on dentin removes the mineral substance greatly reducing the adhesion and must therefore be absolutely avoided. Carefully remove temporary restoration and temporary cement dissolving it with pure acetone or ethanol. It is advisable to use one solution composed of Clorexidina 0,2% + Cetrimide 0,2% for disinfection.

Check the proper fitting of the restoration and the interproximal contacts. If a silicone fit-checker material has been used, carefully remove the residues with specific products, otherwise, adhesion can be impaired. In case of glass-ceramic restorations, adjust the occlusion only after the cementation procedure is accomplished to avoid possible fractures.

Pre-treatment of metal restorations

Follow the indications of the manufacturer of the restoration

materials. If no indications are specified, the following procedure is suggested:

1. Air-abrade the metal surface to be cemented with aluminum oxide 50-100 µm in diameter;
2. Put the prosthesis in pure ethanol or distilled water in an ultrasonic bath for 2 minutes;
3. Wash with water spray and dry with clean air.

Pre-treatment of lithium disilicate and other glass-ceramic restorations

Follow the indications of the manufacturer of the restorative materials. If not specified, the following procedure is suggested:

1. Use hydrofluoric acid 9% for 30-60 seconds to etch the glass-ceramic surfaces;
2. Rinse thoroughly with water spray for 20 seconds and dry with a jet of air free of water and oil (residual surface acidity can be neutralized with sodium bicarbonate 5% solution).
3. It is advisable the use of a proper silane.

Pre-treatment of zirconia and other poly-crystal ceramic restorations

If no indication from the manufacturer of the restorative material are available, the following procedures are suggested:

Procedure 1

1. Air-abrade the restoration surface with aluminum oxide 50-80 µm in diameter at 3,5 bar and to a distance of 10 mm.
2. Immerse the restoration in distilled water or pure ethanol in an ultrasonic bath for 2 minutes in order to remove the residual aluminum oxide particles.
3. Wash and dry the restoration with a clean air stream.

Procedure 2

1. With a sandblaster perform tribo-chemical treatment of the ceramic surface using for example CoJet or Rocatec sand, be sure to direct the air stream orthogonally to the surface.

2. Blow-off the residues of the blasting material with a jet of air free of water and oil. 3. If needed, apply a suitable silane or primer on the treated surface accordingly to the instructions of the manufacturer of these materials.

Pre-treatment of composite restorations

The composite or hybrid ceramic restorations should be treated following the manufacturer's indications. If not specified, the following procedure is suggested:

1. Air-abrade the restoration surface to be cemented with aluminum oxide 50-100 µm in diameter at 2 bar of pressure.
2. Immerse the restoration in distilled water or pure ethanol in an ultrasonic bath;
3. Dry with a jet of air free of water and oil.

Pre-treatment of root canal posts and dowels

Endodontic posts and dowels should be treated following the manufacturer's indications. If not specified, the following procedure is suggested:

1. Clean the fiber post by applying 32-37% orthophosphoric acid gel for 60 seconds or by immersion in ethanol for 5 minutes.
2. Wash with water spray and dry with a jet of air free of water and oil.

If Overfibers fiber posts are used, it is not necessary to apply a silane, but it is advisable to etch them with a phosphoric acid gel 32 - 37% for 60 seconds, then wash with abundant water spray and air-dry. Hi-Rem+ and Over Post+ have the AAS (Active Adhesion System) that increase the adhesion with OverCEM SA luting cement after phosphoric acid-etching treatment.

In case of metal posts or dowels, an air-abrasion treatment of the post surface with aluminum oxide 50-100 µm in diameter followed by water spray cleaning is recommended.

Pre-treatment of screwed implant prosthodontic abutments

OverCEM SA Opaque is intended for screwed-type implant prostheses when a strong connection is required between the implant abutment and the prosthesis material (usually zirconia, glass-ceramic, or composite).

Even if sandblasting of the implant abutment components may improve adhesion, normally a steam-cleaning or pure ethanol cleaning followed by an accurate air-drying is enough to obtain a strong bonding: any treatment that exceeds the normal surface cleaning is optional and must be evaluated in accordance with the instructions of the manufacturer of the implant component.

However, OverCEM SA is not intended for cemented-type implant prosthodontics!

In this case, do not use OverCEM SA because, if needed, excessive adhesion would not allow the detachment of the prosthesis from the implant.

Pre-treatment of metal Maryland and inlay/onlay bridges

Air-abrasion aluminum oxide 50-100 µm at the pressure of 2 bar at the distance of 5-10 mm should be performed on metal surfaces. It is advisable to place the prosthesis in ethanol then in ultrasonic bath with distilled water for at least 5 minutes before cementation in order to eliminate the residual aluminum oxide particles. Otherwise, remove the residues of the sandblasting material with a jet of air free of water and oil. The tribo chemical treatment with CoJet or Rocatec sand followed by a silane or primer application is optional and subjected to the clinician judgement.

Polymerization times

The working and setting times are influenced by the ambient and oral temperature. As for the composite cement, the hardening of OverCEM SA slows down considerably at room temperature and accelerates in the oral cavity. OverCEM SA is a dual-cured material and is therefore sensitive to daylight and artificial light. The working time is reduced in case of handling

under a white light operating lamp: reduce the light intensity to avoid unwanted polymerization. (Tab. 2)

NOTE:

times indicated above are referring to initial temperature of cement 23°C/73.4F.

If cement is stored in a cooler, it is recommended to take it out at least one hour before its use.

Polymerization of OverCEM SA: general recommendations

Two polymerization options are available for OverCEM SA cements: light- and self-curing. OverCEM SA possess a very strong photo-polymerization capability.

In case of translucent restorative materials (composites, ceramics, hybrid ceramics, resins etc.) it is advisable to light-cure the cement through the restoration. For translucent posts, polymerize directly on the post head once it has been properly shortened.

It is advisable to use a dual-band polymerization lamp with an emission ranging from 395 to 480 nm. A minimum power density of 1000 mW/cm² is required to polymerize OverCEM SA in normal conditions. Some lamps have a booster emission option (normally >2000mW/cm²) that increases the light output power. For thicker restorations (>1mm) or endodontic posts the use of the booster-type emission is recommended.

In case of opaque restorative materials (metal, carbon, porcelain-fused-to-metal, etc.) the cement will polymerize mainly by self-curing reactions. Wait for the declared setting time before starting any other procedure that could generate movements or stresses on the restoration.

With opaque restoration, it is advisable to perform light-curing on their margins or where it is likely that even a small amount of curing light can penetrate.

Use of the double-syringe

1. Select a new syringe of OverCEM SA according to the clinical need and prepare it for the use.

2. At the first application remove the syringe cap and dispose it in accord to the local regulations.
3. Before inserting the new T-mixer tip, bleed the syringe so that both base and catalyst pastes are trimmed at equal amounts.

Warning! Make sure that base and catalyst pastes do not come in contact near the exit holes during this operation!

4. Do not re-use the blue closing cap to seal the syringe.
5. Select the T-Mixer tip accordingly to the clinical needs. Verify the wedge-shaped protrusion on the blue ferrule located on the T-mixer ring and insert it into the corresponding triangular notch of the syringe. Then rotate the ring 90°. The T-Mixer Regular with conical tip is used to apply the cement in easily accessible areas, whereas to inject OverCEM SA into the root canal, use the specific T-Mixer Wide + Endo tips.

Cementation of indirect restorations

The tooth surfaces to be cemented should be clean but not excessively dried. Avoid contamination with water, blood, saliva and crevicular fluid before cementation and during its hardening phase.

1. Evenly distribute OverCEM SA on the restoration surfaces by pressing the pistons and extruding the material in slight excess. Keep the T-Mixer tip immersed in the material during the entire application to prevent the inclusion of air bubbles.
 2. Apply the restoration on the substrate and stabilize it with a constant finger or medium strength pressure of patient bite.
 3. Wait at least 60 seconds before attempting any light curing: this will allow the chemical reaction with the dentin and enamel;
- Light-cure the cement excess with a short burst of 1-2 seconds on each tooth side to make it solid enough to be easily removed. Use a dental floss to carefully clean the interproximal or connector areas.

4. Continue the light-curing procedure at the specified times (Tab. 2).
5. Leave the used T-Mixer on the syringe until the next application: it will replace the cap until the next use avoiding the contamination between pastes. Store the syringe in a place protected from light and heat.
6. Finish and polish the restoration with the appropriate rotating instruments (rubber cups and wheels, brushes, etc.). Check the occlusion to avoid any interference with the opposing teeth.

Fiber post and dowels cementation

1. Mount the Endo tip in the T-Mixer Wide and insert it on the double syringe.
2. Introduce the Endo tip deeply into the root canal then start to inject OverCEM SA.
3. Keep the tip immersed in the cement and slowly move the tip upwards in accordance with the canal filling speed. Do not extract the tip from the cement until the root canal is full. The cement injection must be properly carried out to minimize the inclusion of air bubbles.
4. Translucent fiber posts: if present, remove the colored post head using a diamond bur with cooling spray, then wet the coronal surface with a small quantity of cement: this will allow a better light passage in the subsequent curing phases;
5. Insert the post in the canal already filled with OverCEM SA and gently push it to the bottom by applying a rotating movement for a better seating. Using a brush or a sponge, evenly distribute the cement on the tooth and outer post surface, in order to create a bonding layer for the core build-up composites that will be subsequently applied. Remove all the cement excesses.
6. Light-polymerize for at least 60 seconds making sure the light flow is parallel to the post axis. If available, use the booster light emission for at least 40 seconds.

WARNING! During polymerization, avoid any contact between the lamp and the post head it may impede correct adhesion!

Final checks

Use optical magnification systems to search and eliminate residuals of cement around the teeth and inside the gingival sulcus.

Proceed to radiographic examination if any doubt persists.

Hygiene and disinfection

Clean OverCEM SA syringes after the use with detergents and surface disinfectants commonly used in the dental practice.

Precautions

Avoid exposure of the cement to direct bright light or humidity during cement application to avoid setting time shortening.

Conservation and expiry

OverCEM SA should be stored between 2°C/36°F-25°C/77°F. After the first use, it is recommended the storage from 2°C/36°F to 8°C/46°F for a better conservation of OverCEM SA. For better results do not use the OverCEM SA cements after their expiry date.

Information about the product

Overfibers is the only source of correct information for using OverCEM SA.

It is recommended to keep the instructions for use during the entire duration of use of the product.

Warranty

Overfibers warrants this product to be free from defects in materials and workmanship. Overfibers do not offer further warranties, including any implied warranties or of merchantability or fitness for particular purposes. The user is responsible for determining the suitability of the product in the individual applications. If this product is defective within the warranty period and Overfibers ascertains an anomaly of the same, it will be replaced.

Limitations of responsibility

Except where otherwise indicated by law, Overfibers will not be held responsible for any loss or damage resulting from this product, direct or indirect, special, incidental or consequential, whatever the established theory, including warranty, contract, negligence or direct responsibility.

Implementation of the instructions for use: 06/2020, Rev. 02

CE
0425

Tab. 2

OverCEM SA	Version	Working Time	Setting Time	Working Time	Setting Time	Light-curing Occlusal surface	Light-curing any other surface	Light-curing of posts
		35°C/95F		23°C/73.4F				
	Universal	01:30	04:00	03:00	06:00	00:20	00:20	00:40
	Translucent	01:30	04:00	03:00	06:00	00:20	00:20	00:40
	Opaque	01:30	04:00	03:00	06:00	00:30	00:30	00:60

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